
English	INSTRUCTIONS FOR USE: Southern Implants® Instrument Tray
Español	INSTRUCCIONES DE USO: Bandeja para instrumental Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Vassoio per Strumenti
Français	MODE D'EMPLOI : Plateaux à instruments de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Instrumententablett
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® tabuleiros de instrumentos
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Δίσκος οργάνων
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Instrumentbricka
Brazilian	INSTRUÇÕES DE USO: Southern Implants® Bandeja de instrumentos
Hrvatski	UPUTE ZA UPOTREBU: Tacna za instrumente Southern Implants®
Dansk	BRUGSANVISNING: Southern Implants® Instrumentbakke
Suomi	KÄYTTÖOHJEET: Southern Implants® uudelleenkäytettävät instrumentit
Nederlands (NL)	GEBRUIKSAANWIJZING: Southern Implants® instrumententray
Magyar	HASZNÁLATI ÚTMUTATÓ: Southern Implants® műszertálca
Lietuviškai	NAUDOJIMO NURODYMAI: „Southern“ implantai® Instrumentų dėklas
Polski	INSTRUKCJA UŻYTKOWANIA: Southern Implants® Taca na instrumenty
Română	INSTRUCȚIUNI DE UTILIZARE: Tavă instrumentar Southern Implants®
Slovenščina	NAVODILA ZA UPORABO: Pladenj za instrumente Southern Implants®
Українська	ІНСТРУКЦІЇ З ЕКСПЛУАТАЦІЇ. Лоток для інструментів Southern Implants®
Русский	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ: Лоток для инструментов Southern Implants®
Български	ИНСТРУКЦИИ ЗА УПОТРЕБА: Тави за инструменти Southern Implants®
Česky	NÁVOD K POUŽITÍ: Podnosy na nástroje Southern Implants®
Eesti	KASUTUSJUHEND: Southern Implants®-i instrumendialus
Íslenska	NOTKUNARLEIÐBEININGAR: Southern Implants® áhaldabakki
Latviski	LIETOŠANAS INSTRUKCIJA: Southern Implants® instrumentu paplāte
Norsk	BRUKSANVISNINGER: Southern Implants® Instrumentbrett
Slovenčina	NÁVOD NA POUŽITIE: Prístrojová podložka Southern Implants®
Türkçe	KULLANMA TALİMATI: Southern Implants® Alet Tepsisi
Српски	УПУТСТВО ЗА УПОТРЕБУ: Southern Implants® подметач за инструмент

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customer care@southernimplants.com

Description

The device 'instrument tray' is a reusable rigid container or organising tray intended for use in health care facilities for the purpose of organising, containing, and transporting reusable medical devices for sterilisation. It is composed of multiple pieces, designed to be integrated into a single unit which contains and protects the interior components during steam sterilisation. Each tray consists of two components: a base tray and a lid, with some trays containing a third component - an internal individualised insert tray. The tray base and lid are perforated for steam sterilisation, with the insert tray allowing steam to penetrate through the retention sockets (grommets) which are made from silicone. The internal insert tray and base has the ability to hold individualised pieces and accessories which include dental tools, drills and ratchets/wrenches.

The lid is made from Radel R-5800 and is tinted a transparent blue or grey. The base and insert tray are made of Radel R-5000 and is white/light grey in colour. The Radel material is a polymer resin. The base and insert trays are layered with medical grade silicon layers to locate and hold the various tools and instruments.

Intended use

The instrument trays are used to hold dental tools and instruments in place during transport, steam sterilisation and storage. The instrument trays are reusable devices. The instrument trays are medical devices.

Indications for use

The Southern Implants' Instrument Tray is designed to hold various dental surgical drills and tools in order to organise, steam sterilise, and transport the instruments between uses. The tray is to be enclosed in an FDA cleared steam sterilisable wrap and sterilised in an FDA Cleared steriliser for one of the following cycles:

- Pre-vacuum steam – At 132°C for 4 minutes with a 20 minutes dry time.
- Pre-vacuum steam – At 135°C for 3 minutes with a 20 minutes dry time.

The trays are not intended for sterilisation of non-porous loads. The trays are recommended not to be stacked during sterilisation. The Complete Surgical Trays represent the worst case validated load due to number of components (Large: 25 instruments; Medium: 90 instruments; Small: 47 instruments) and the weight (Large: 752 grams; Medium: 672 grams; Small: 339 grams). Southern Implants (Pty) Ltd does not make any lumen claims for the Southern Implants' Instrument Trays.

Size (L x W x H)	Product Code	Number of Instruments	Weight of Full Tray (g)
Large 26.8 x 14.7 x 5.5 cm	CH-I-ZYG	19	752
	I-ZYG-1	25	746
Medium 18.7 x 13.5 x 5.5 cm	I-HEX-EG	90	672
	I-IT-EG	49	541
	I-TRI-NEX-EG	60	547
	I-DC-EG	53	540
	I-INT-HEX-EG	46	555
	I-MAX-EG	50	546
Small 14.8 x 9.5x 5.5 cm	I-IV-EG	47	339
	I-SP-EG	44	310
	I-SP-GS1-EG	47	339
	I-SP-GS2-EG	47	339
	I-PROS-EG	25	280
	I-PROS-MINI	7	147

Intended user

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and experienced implant users.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theater or a dentist consultation room.

Intended patient population

This device is used in the dental restoration of partially or fully edentulous patients in the upper or lower jaw. Restorations may comprise of single teeth, partial or full bridges, and may be fixed or removable.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants' range there are 11 implant systems. The implant code and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

Table A – Corresponding instrument tray with implant system

Product Code	Implant System
CH-I-ZYG	Straumann Zygomatic Implant System
I-ZYG-1	Southern Implants Zygomatic Implant System
I-HEX-EG	External Hex and PROVATA® Implant Systems
I-IT-EG	Internal Octagon (IT) Implant System
I-TRI-NEX-EG	TRI-NEX® Implant System
I-DC-EG	Deep Conical Implant System
I-INT-HEX-EG	Internal (M-Series) Hex Implant System
I-MAX-EG	MAX Implants including TRI-MAX®, MAX, PROMAX®, MAXIT®
I-IV-EG	INVERTA® Deep Conical & INVERTA® External Hex Implants
I-SP-EG	Single Platform (SP1) Implant System
I-SP-GS1-EG	Single Platform (SP1) Implant System
I-SP-GS2-EG	Single Platform (SP1) Implant System
I-PROS-EG	Applicable to all Southern Implant systems
I-PROS-MINI	Applicable to all Southern Implant systems

Clinical benefits

Through this procedure patients can expect to have their missing teeth replaced and/or crowns restored.

Storage, cleaning and sterilisation

Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device, it should be thoroughly inspected and tested to determine its suitability for reuse.

NOTE: during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

NOTE: always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermodisinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

Inspection

Do a visual inspection of the items to check for damages on instruments.

Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the restoration prior to use:

1. prevacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

Contraindications

The contraindications for this system are not dissimilar to those of dental implant therapy. Contra-indications to implant therapy include:

- patients medically unfit for oral surgical procedures.
- where an inadequate number of implants can be placed limiting the functional support of the prosthesis.
- patients under the age of 18 years.
- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients with drug or alcohol abuse dependencies.
- patients undergoing chronic high dose steroid therapy.
- patients undergoing anti-coagulant therapy.
- patients with metabolic bone disease.

- patients undergoing radiotherapy treatment.
- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, or iridium allergies or hypersensitivities.

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Products must be secured against aspiration when handled intraorally. Aspiration of products may lead to infection or unplanned physical injury.

Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone. Implant failure increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue, leading to diminished healing capacity.

It is important to be aware and avoid damage to vital structures such as nerves, veins and arteries. Injuries to vital anatomical structures may cause serious complications such as injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve. Failing to identify actual measurements relative to the radiographic data could lead to complications.

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had orofacial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimising the trauma to the host tissue increases the potential for successful osseointegration.
- electrosurgery should not be attempted around metal implants as they are conductive.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information for the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The side effects of the use of the system are not dissimilar to those of dental implant therapy. Possible side effects to implant therapy include:

- pain
- swelling
- phonetic difficulties
- gingival inflammation

Less common but more persistent symptoms include, but are not limited to:

- allergic reaction(s) to implant and/or abutment material
- breakage of the implant and/or abutment

- loosening of the abutment screw and/or retaining screw
- infection requiring revision of the dental implant
- nerve damage resulting in permanent weakness, numbness, or pain
- histologic responses with possible macrophage and/or fibroblast involvement
- fat emboli formation
- loosening of the implant requiring revision surgery
- perforation of the maxillary sinus
- perforation of the labial and lingual plates
- bone loss possibly resulting in revision or removal of the implant.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK “peel-back” lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com.

Materials

Material type Radel R-5800 (lid), Radel R-5000 (base and insert),
Medical Grade Silicon (Grommet inserts)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.




Basic UDI

Product	Basic-UDI Number
Basic-UDI for Instrument Trays	6009544039118P

Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue
 CAT-2005 - IT Implants Product Catalogue
 CAT-2010 - Osseointegrated Fixtures Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue
 CAT-2093 - Single Platform Implants Product Catalogue

Symbols and warnings

 Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	 CE mark 2797	 Prescription device*	 Sterilised using irradiation	 Non-sterile	 Use by date (mm-yy)	 Do not reuse	 Do not resterilise	 Catalogue number	 Batch code	 Medical device	 Authorised representative in the European Community
 Authorised representative for Switzerland	 Date of manufacture	 Magnetic Resonance conditional	 Magnetic Resonance safe	 Single sterile barrier system with protective packaging inside	 Single sterile barrier system	 Consult instruction for use	 Caution	 Keep away from sunlight	 Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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